

# CL10

## MAGNETIC COVER LIFTER

unique on the market usable in 4 different ways!



Universal model with rectangular base ideal for lifting, moving and repositioning of the majority of covers and grilles in cast iron and ferrous materials, the handle is telescopic and folding to allow ergonomics and minimal space for transport and stocking. The exclusive configuration gives versatility of use in four ways to fit the majority of lifting/repositioning operations. The plate with 16 powerful neodinium magnets assures a lifting force of 300 Kg with a weight of only Kg 7. **The use is simple, fast and intuitive.**

### CHARACTERISTICS

WEIGHT	Kg 7	BASE DIMENSION	cm 21 x 16
LENGTH IN EXTENSION			cm 70
LENGTH RETRACTED/FOLDED			cm 55 x 16



• **ANTI MAGNETIC PROTECTIVE SACK** Code A221006  
Accessory of the magnetic cover lifters, prevents unwanted contacts between the magnetic plate with ferrous objects that could be difficult to separate and render normally manageable the tools. (data sheet at pag. 11)

1



#### MANUAL LIFTING WITH MAGNETIC PLATE

CL10 is placed on the cover in the position nearest to the frame giving a solid blockage, the telescopic handle allows to manually lift and move the majority of covers, to detach the tool is enough to give a strong side stroke.

**The operation is simple and fast.**

2




#### LIFTING WITH APS90 - APS80 LEVERS

Lifting using CL10 with the folding levers the efforts and possible accidents on the back are totally eliminated, the operator weight becomes "friend" thanks to the downward push action that stretches the vertebrae and back muscles,

**The operations are simple, fast and safe.**

3



#### LIFTING WITH TWO MAGNETIC PLATES BLOCKED IN PARALLEL AND APS90 - APS80 LEVERS

If cover is too heavy or too blocked by dirt or too grooved with minimal magnetic contact is possible to use two CL10 in parallel united by a small connection "bridge", the contact force is doubled, the operator weight becomes "friend" thanks to the downward push action that stretches the vertebrae and back muscles, the plates will be then detached one at a time.

**The operations are simple, safe and fast.**

4



#### MANUAL LIFTING WITH HANDLE BLOCKED HORIZONTALLY TO LIFT HINGED ROUND COVERS

By simply changing the position of a pin is possible to block in horizontal position the telescopic handle. In this configuration the magnetic plate becomes a practical rigid long handle to lift and reverse the covers that are hinged to the frame.

**The operation is simple, fast and safe.**